TEST DE ANTÍGENOS AUTODIAGNÓSTICO COVID-19

CE 0123

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Product Service

EC Certificate

EC Design-Examination Certificate

Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6)
(Devices for self-testing)

No. V9 061317 0006 Rev. 00

Manufacturer: **Xiamen Boson Biotech Co., Ltd.**
90-94 Tianfeng Road
Jimei North Industrial Park
361021 Xiamen, Fujian
PEOPLE'S REPUBLIC OF CHINA

Product: **In Vitro diagnostic devices for self testing**

Model(s): **Rapid SARS-CoV-2 Antigen Test Card**

Parameters:

Model Name:	Model No.:
Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-2
Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-4
Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-6

The Certification Body of TÜV SÜD Product Service GmbH declares that a design examination has been carried out on the respective devices in accordance with IVDD Annex III (6). The design of the devices conforms to the requirements of this Directive. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:V9_061317_0006_Rev_00

Report No.: 713210321

Valid from: 2021-04-01

Valid until: 2022-05-26

Date, 2021-04-01

Christoph Dicks
Head of Certification/Notified Body

Page 1 of 1

TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany

TÜV®

Declaration of Conformity

Manufacturer	Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P. R. China.
European Representative	Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Product	Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-2 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-4 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-6
Classification	Self-Test IVD
Confirmaty Assessment Route	Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III section 6

We herewith declare that above mentioned products meet the provisions of the council Directive 98/79/EC for medical devices. All Supporting Documentation is Retained under the premises of the manufacturer. We have sole responsibility for issuing the Declaration of Conformity.

Standard Applied	ISO 13485:2016	EN ISO 14971:2012
	EN ISO 18113-1:2011	EN ISO 18113-4:2011
	ISO 15223-1:2016	EN 62366-1:2015
	EN 13612:2002	EN ISO 23640:2015
Notified Body	TÜV SÜD Product Service GmbH, Ridlerstraße 65, 80339 Munich, Germany. Identification number: 0123	
(EC) Certificate(s)	No. V9 061317 0006 Rev. 00	
Start of CE-Marking	2021-04-01	
Place, Date of Issue	Xiamen, 2021-04-01	

Signature Changgong Zhang
(Signed By Boson Representative)
Name: Changgong Zhang
Title: General Manager

CE
0123

Certificate

No. Q5 061317 0005 Rev. 00

Holder of Certificate: **Xiamen Boson Biotech Co., Ltd.**
 90-94 Tianfeng Road
 Jimei North Industrial Park
 361021 Xiamen, Fujian
 PEOPLE'S REPUBLIC OF CHINA

Facility(ies): Xiamen Boson Biotech Co., Ltd.
 90-94 Tianfeng Road, Jimei North Industrial Park, 361021
 Xiamen, Fujian, PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate: **Design and Development, Production and Distribution of In Vitro Diagnostic kits for detection of Infection Diseases, Tumour Markers, Drug Abuse, Hormones, Cardiac Markers and Related Biomaterial**

Applied Standard(s): EN ISO 13485:2016
 Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
 DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: SH1807513

Valid from: 2018-10-31

Valid until: 2021-10-30

Date, 2018-09-27

Stefan Preiß

TÜV SÜD
 ZERTIFIKAT ◆ CERTIFICATE ◆ 認 證 證 書 ◆ CERTIFICADO ◆ CERTIFICAT

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English RAPID SARS-COV-2 ANTIGEN TEST CARD INSTRUCTION GUIDE FOR ANTERIOR NASAL SWAB SPECIMENS For Use at Home

REF	1N40C5-2	For 1 Test/Box
REF	1N40C5-4	For 5 Tests/Box
REF	1N40C5-6	For 20 Tests/Box

Please follow the instruction label carefully.
INTENDED USE:
 Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based on one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 antigen in swab specimens from individuals suspected of COVID-19 within the first seven days of symptom onset. Rapid SARS-CoV-2 Antigen Test Card shall not be used as a basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be tested in a clinical laboratory setting.

SUMMARY:
 The novel coronavirus binds to the β-glycan. COVID-19 is an acute respiratory infectious disease. People are generally asymptomatic. Currently, the patients infected by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Loss of taste, loss of smell, sore throat, runny nose, and conjunctivitis are also common.

MATERIALS PROVIDED:

Components	For 1 Test/Box	For 5 Tests/Box	For 20 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pack)	1	5	20
Swab	1	5	20
Control swab	1	5	20
Instructions for use (this booklet)	1	1	1

PERFORMANCE (SENSITIVITY AND SPECIFICITY)
 Rapid SARS-CoV-2 Antigen Test Card was compared to the confirmed clinical diagnosis. The study involved 156 samples.
Sensitivity: 99.23%
Specificity: 99.23%

INTERPRETATION:
 A visibly dark spot demonstrated that:
 - 99.23% of individuals carried out the test without requiring assistance.
 - 97.87% of the different types of results were interpreted correctly.

LIMITATIONS:
 None of the following substances at the tested concentration showed any interference with the test:
 White Blood Cell: Alkaloid 1%, Mucin 2%,
 Phosphate: 10%, Oxymetazolin 0.004%,
 Menthyl 0.15%, Camphor 1%, Benzocaine 0.15%,
 Fluorocaine Phosphate 0.1%, Eucalyptol 0.1%,
 Osmolality: 300 mOsm/kg, Sodium chloride 0.5%, Human Antinuclear Antibody (HAMA): 0.01ng/ml.

IMPORTANT INFORMATION BEFORE THE EXECUTION:
 1. Read the instruction guide carefully.
 2. Do not use the product beyond the expiration date.
 3. Do not use the product if the pouch is damaged or the seal is broken.
 4. Store the test device at 15 to 30°C in the original sealed pouch. Do not freeze.
 5. The test device should be used at room temperature (15°C to 30°C). If the product has been stored in a cool area (less than 15°C), leave it at normal room temperature for 30 minutes before using.
 6. Handle all specimens as potentially infectious.
 7. Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
 8. Use the swabs included in the test kit to ensure optimal performance of the test specimen.
 9. Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretions) with the swab, especially for anterior nasal sampling.
 10. Blow the nose several times before collecting the specimen.
 11. The specimens should be tested as soon as possible after collection.
 12. Apply the excess of test specimen to the specimen tube.
 13. Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
 14. When used as intended, there should not be any contact with the extraction buffer, a case of contact with the eyes, mouth or other parts, rinse with clear water if an irritation persists, consult a medical professional.
 15. Children under 14 years of age should be assisted by an adult.

LIMITATIONS:
 This test is to be used exclusively for the qualitative detection of SARS-CoV-2 antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as a part of this test.
 1. Proper specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretions) with the swab, especially for anterior nasal sampling.
 2. Apply the excess of test specimen to the specimen tube.
 3. Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
 4. When used as intended, there should not be any contact with the extraction buffer, a case of contact with the eyes, mouth or other parts, rinse with clear water if an irritation persists, consult a medical professional.
 5. Children under 14 years of age should be assisted by an adult.

PREPARATION:
 1. Check the test kit contents. Make sure that nothing is damaged or broken.
 2. Tear off the seal.
 3. Blow the nose several times before collecting specimen.
 4. Wash hands.

DISPOSAL:
 The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

PROCEDURE:
 This test is suitable for people of all ages. The recommended operation are again from 14 to 90. Children under 14 years of age should be tested by an adult. Do not continue the test if the child tests any pain.
 Rotate the lid of your extraction buffer bottle and be careful not to spill any of the liquid.
 Open! Open! Open! Open! Open! Open!

1. Squeeze the extraction buffer out of the bottle into the extraction tube.
CAUTION: Avoid touching the nozzle against the tube.

2. Find the swab in the sealed wrapper in front of you. Identify the soft, fabric tip of the swab.

3. Peel open the swab packaging and gently take out the swab.
CAUTION: Never touch the soft, fabric tip of the swab with your hands.

4. Carefully insert swab into one nostril.
CAUTION: Results after 15-20 minutes may not be accurate. The used device will be disposed of with normal household waste in accordance with the applicable local regulations.

5. Pinch extraction tube with fingers and remove the solution from swab as much as possible.

6. Install the nasette cap on the sample extraction tube tightly.

7. Blow the kit components to room temperature before testing. Open the pouch and remove the card. Bring the card to a flat and level surface.
CAUTION: Once opened, the test card must be used immediately.

8. Invert the test card and add 3 drops (15 µl) of test specimen into the specimen tube.
CAUTION: The extraction of air bubbles in the specimen well (S) must be avoided.

9. Read the results at 15-20 minutes.
CAUTION: Results after 20 minutes may not be accurate. The used device will be disposed of with normal household waste in accordance with the applicable local regulations.

INTERPRETATION OF RESULTS:
QUALITY CONTROL:
 The control line is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagent is not expired.
FREQUENTLY ASKED QUESTIONS (FAQ):
 1. How does the detection work?
 2. When should I test myself?
 3. What can affect my test result?
 4. How should I store the test kit?
 5. How can I dispose of the product?

CE Marking:
 The CE mark is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagent is not expired.

WARNING:
 H317: Warning! Liquid component may cause an allergic skin reaction.

Manufacturer: Xiamen Bosen Biotech Co., Ltd.
 90-94 Taining Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P.R. China.
Authorized Representative: Lohus N.V.
 Koningrijks Jachtvliet 10, 1e verd. 2595AA, The Hague, Netherlands.

Version 5.1 Date: April 29th, 2021

Accessories	Manufacturer	EC-Representative	CE Mark
Swab A	Jiangsu Changfeng Medical Industry Co., Ltd. Tongzhu Town, Guangling District Yangzhou 225109 P.R. China.	Lohus Service & Consulting GmbH Frieden-Straße 30/32/24 Hörsburg Germany	enlist 3/REG/EEC
Swab B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongchang Street Jizhou District Dalian 116100 Liaoning, China.	CAC Medical Devices & Drugs S.L. C/ Horacio Lango No.18, CP 29006, Málaga, Spain	enlist 3/REG/EEC
Swab C	Zhejiang Gongchang Medical Technology Co., Ltd. No. 809, Xueyuan Road, Hangzhou 311802 Taishou, Zhejiang, P.R. China	Shanghai International Holding Corp. (Germany) Frieden-Straße 30/32/24 Hörsburg Germany	enlist 3/REG/EEC
Swab D	Jiangsu Hanheng Medical Technology Co., Ltd. 16-28/11 North Guangyuan Road, Taijiang District, 213017, Changzhou, Jiangsu, China.	Lohus Lebenservice GmbH Kochstr. 4/7877, Wilsch, Germany	enlist 3/REG/EEC

EXPLANATION FOR SYMBOLS:
 IVD: In vitro Diagnostic Use
 LOT: Lot Number
 REF: Reference
 CE Mark: CE Marking
 H317: Warning! Liquid component may cause an allergic skin reaction.

Accessories	Manufacturer	EC-Representative	CE Mark
Swab A	Jiangsu Changfeng Medical Industry Co., Ltd. Tongzhu Town, Guangling District Yangzhou 225109 P.R. China.	Lohus Service & Consulting GmbH Frieden-Straße 30/32/24 Hörsburg Germany	enlist 3/REG/EEC
Swab B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongchang Street Jizhou District Dalian 116100 Liaoning, China.	CAC Medical Devices & Drugs S.L. C/ Horacio Lango No.18, CP 29006, Málaga, Spain	enlist 3/REG/EEC
Swab C	Zhejiang Gongchang Medical Technology Co., Ltd. No. 809, Xueyuan Road, Hangzhou 311802 Taishou, Zhejiang, P.R. China	Shanghai International Holding Corp. (Germany) Frieden-Straße 30/32/24 Hörsburg Germany	enlist 3/REG/EEC
Swab D	Jiangsu Hanheng Medical Technology Co., Ltd. 16-28/11 North Guangyuan Road, Taijiang District, 213017, Changzhou, Jiangsu, China.	Lohus Lebenservice GmbH Kochstr. 4/7877, Wilsch, Germany	enlist 3/REG/EEC

Manufacturer: Xiamen Bosen Biotech Co., Ltd.
 90-94 Taining Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P.R. China.
Authorized Representative: Lohus N.V.
 Koningrijks Jachtvliet 10, 1e verd. 2595AA, The Hague, Netherlands.

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Letter of Authorization

To whom it may concern:

We **Xiamen Boson Biotech Co., Ltd.**, a manufacturer of various types of Diagnostic Kits, having factories at 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, P.R.C. do hereby confirm that **PACIFIC MEDICAL CARE S.L.**, address: Doctor Fleming 31,10 planta, 28036 Madrid, Spain, is our authorized local agent/distributor to import, register, commercialize and distribute our following product in Spain. This authorization will automatically expire when Boson authorizes an exclusive distributor in Spain.

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